

i.Need4

Your eyes need I.NEED4



EMPOWER YOUR EYELID WITH 4MHZ RF SYSTEM

- Prompt minor wrinkle removing effect especially around the eyes.
- Prevent for eyelid sagging problem if you continue this treatment, no needs blepharoplasty surgery.
- Lower the eye strain from computers, televisions and other devices you use daily.



I.NEED4 STRENGTH

- 4MHz high frequency generates diathermy heat in a short time
- Electrode is designed for safe and convenient treatment for highly sensitive small area like upper and lower eyelids, crow's feet
- There is no shot limit, Handpiece can be reusable ► **Economical**
- Powerful Synergy effect by treating 10 to 20 min. treatment after other full face care devices
- Prompt elasticity increase, wrinkle reducing effect
- Treatment area: focusing on around eyes, perioral, sagging pore and lips

MECHANISM

4.0MHz Monopolar RF treatment provides prompt wrinkle removing effect by generating neo-collagen and relieving tension and stiff muscles around the eyes.



Old and disarranged collagen
in dermis make wrinkles on the skin

4MHz High Frequency RF very extremely
quickly increase heat to 39°C to 45 °C

1st Prompt skin firm effect
2nd New Collagen and Elastin are generated


I NEED 4 PREVENTS LID PTOSIS BY AGING

I needs 4 gives more precise lifting effect with its specially designed electrode and 4MHz high frequency. When the electrode contacts with skin, it flows high frequency (4MHz) into skin and increase temperature in dermis and it restructures the collagen and change to compact.

Lid Ptosis is most common in older adults because of the natural aging process.

The levator muscle is responsible for lifting the eyelid. As you age, that muscle can stretch and, as a result, cause the eyelid to fall. Keep in mind, though, that people of all ages can be affected by this condition.

In some cases, this can affect your vision. Since ptosis can affect your vision, you should take it seriously.



Distance between upper lid margin and light reflex (MRD)



Mild Ptosis (2mm of droop)



Moderate Ptosis (3mm of droop)





Severe Ptosis (4mm or more droop)

I NEED 4 IMPROVES EYE HEALTH

Prolonged viewing of digital devices such as smartphones, computer screen, TV and etc. can be cause of tired eye and aggravate eye health. I need 4 relieves tension and stiff muscles focusing around the eyes with RF diathermy heating. During I need 4 care oxygen level increase in eye so it makes patient feel fresh and helpful for improving eye health.

Eye edema, dark circle, scheroma by accumulation of fatigue and tired eyes by eye muscle tension and hypofunction of lacrimal glands can be cured by the eye-care



Technical Specification

Name of Brand : INEED4	Model Number : INE-7200
Name of Product : RF Diathermy device	Rated Voltage : 100-240Vac , 50/60Hz
Size of the main body : 310*231*142mm / 7" Touch Screen	Weight : 3kg

